

Simple Stuff Works



SLEEPING & POSITIONING SYSTEM

Distributed by
medifab

SCRIPTING AND ORDERING GUIDE

Client details

Client/User name:

DOB:

Date:

Therapist:

Notes

Further Information

CHOOSE A **KIT** TO USE ON **FOAM MATTRESS**

Kit Suggestion

Baby Kit <input type="radio"/>	Child Kit <input type="radio"/>	Teen/Adult Kit <input type="radio"/>
0-5 Years	5-12 Years	12+ Years
<p>Package includes:</p> <ul style="list-style-type: none">• Cot Stabilising Mesh Pack• Cot Temperate Fibre Topper• Fitted Cotton Sheets – Cot Size• Small Laterals• Small Wedges• Small Horseshoe pillow• Simple Stabiliser Size 1• Travel Bag	<p>Package includes:</p> <ul style="list-style-type: none">• Single Bed Stabilising Mesh Pack• Single Bed Temperate Fibre Topper• Fitted Cotton Sheets – Single Bed Size• Medium Lateral Supports• Small Wedges• Small Lateral Supports• Simple Stabiliser Size 2• Small Terry Towelling Pillow• Travel Bag	<p>Package includes:</p> <ul style="list-style-type: none">• Single Bed Stabilising Mesh Pack• Single Bed Temperate Fibre Topper• Fitted Cotton Sheets – Single Bed Size• Medium Lateral Supports• Large Lateral Supports• Medium Wedges• Large Terry Towelling Pillow• Travel Bag
7150-7001-C00	7150-7002-S00	7150-7003-S00

OR REFER TO **PAGE 3** TO **CONFIGURE** YOUR OWN

CONFIGURE YOUR OWN
FOLLOW THE PATHS DEPENDING ON CLIENT'S MATTRESS



Choose your set for foam mattress

Cot Bed	Single Bed	
<p>Set Includes:</p> <ul style="list-style-type: none"> • Stabilising mesh set (includes 2 mesh sheets; 1 bottom, 1 top) • 2 x cot temperate toppers • 2 x cot fitted sheets <p>7150-7021-C00 <input type="radio"/></p>	<p>Set Includes:</p> <ul style="list-style-type: none"> • Stabilising mesh set (includes 2 mesh sheets; 1 bottom, 1 top) • 2 x single temperate toppers • 2 x single fitted sheets <p>7150-7022-S00 <input type="radio"/></p>	

Choose your set for air mattress

Single Bed	
<p>Set Includes:</p> <ul style="list-style-type: none"> • 2 x single temperate toppers • 2 x single fitted sheets <p>7150-7024-S00 <input type="radio"/></p>	




Soft wedges for foam mattress

Soft wedges for air mattress

Wedges with padded covers (pair)		Wedges with padded covers and integrated mesh (pair)	
<p>Small</p> <p>7150-7351-000 <input type="radio"/></p>	<p>Medium</p> <p>7150-7352-000 <input type="radio"/></p>	<p>Small</p> <p>7150-7326-000 <input type="radio"/></p>	<p>Medium</p> <p>7150-7327-000 <input type="radio"/></p>




Lateral supports for foam mattress

Lateral Supports with padded covers (pair)

Small 7150-7201-000	Medium 7150-7202-000	Large 7150-7203-000
		


Lateral supports for air mattress

Lateral Supports with padded covers and incorporated mesh (pair)

Small 7150-7211-000	Medium 7150-7212-000	Large 7150-7213-000
		


Neck support pillow

Promotes neutral neck alignment

<p>(suitable for air & foam mattresses) 300mm x 350mm x 100mm 7150-4050-000</p>	
<p>Spare terry cotton cover 7150-4055-000</p>	


Head Positioning Pillow

Assists with head positioning

Small A: 480mm B: 350mm		Large A: 680mm B: 450mm		
7150-7501-000	<input type="radio"/>	7150-7502-000	<input type="radio"/>	
Spare terry cotton cover 7150-7511-000	<input type="radio"/>	Spare terry cotton cover 7150-7512-000	<input type="radio"/>	


Roll Pillows

Assists with trunk and/or thigh positioning

Small A: 150mm B: 400mm		Large A: 180mm B: 500mm		
With breathable inner & terry cotton outer 7150-7851-000	<input type="radio"/>	With breathable inner & terry cotton outer 7150-7852-000	<input type="radio"/>	
With wipeable incontinence inner & terry cotton outer 7150-7855-000	<input type="radio"/>	With wipeable incontinence inner & terry cotton outer 7150-7856-000	<input type="radio"/>	
Spare terry cotton cover 7150-7857-000	<input type="radio"/>	Spare terry cotton cover 7150-7858-000	<input type="radio"/>	


Horseshoe Pillows

Multi-purpose positioning pillow

Small A: 190mm B: 460mm		Large A: 200mm B: 500mm		
With breathable inner & terry cotton outer 7150-7711-000	<input type="radio"/>	With breathable inner & terry cotton outer 7150-7712-000	<input type="radio"/>	
With wipeable incontinence inner & terry cotton outer 7150-7721-000	<input type="radio"/>	With wipeable incontinence inner & terry cotton outer 7150-7722-000	<input type="radio"/>	
Spare terry cotton cover 7150-7731-000	<input type="radio"/>	Spare terry cotton cover 7150-7732-00	<input type="radio"/>	

Simple Stabiliser

Supports hips and thighs in neutral supine alignment

Size 1	Size 2	Size 3	Size 4
H: 11cm D: 11cm W: 42cm To fit approx leg circumference of 25cm. Suggested age range 2 – 4yrs.	H: 19cm D: 35cm W: 46cm To fit approx leg circumference of 35cm. Suggested age range 2 – 10yrs.	H: 25cm D: 40cm W: 65cm To fit approx leg circumference of 55cm. Suggested age range 5 – 13yrs.	H: 28cm D: 40cm W: 65cm To fit approx leg circumference of 60cm. Suggested age range 13 – adult.
With breathable inner & terry cotton outer 7150-7341-000 <input type="radio"/>	With breathable inner & terry cotton outer 7150-7342-000 <input type="radio"/>	With breathable inner & terry cotton outer 7150-7343-000 <input type="radio"/>	With breathable inner & terry cotton outer 7150-7344-000 <input type="radio"/>
With wipeable incontinence inner & terry cotton outer 7150-7345-000 <input type="radio"/>	With wipeable incontinence inner & terry cotton outer 7150-7346-000 <input type="radio"/>	With wipeable incontinence inner & terry cotton outer 7150-7347-000 <input type="radio"/>	With wipeable incontinence inner & terry cotton outer 7150-7348-000 <input type="radio"/>
			
Spare terry cotton cover 7150-7221-000 <input type="radio"/>	Spare terry cotton cover 7150-7222-000 <input type="radio"/>	Spare terry cotton cover 7150-7223-000 <input type="radio"/>	Spare terry cotton cover 7150-7224-000 <input type="radio"/>
Spare wipeable inner cover 7150-7225-000 <input type="radio"/>	Spare wipeable inner cover 7150-7226-000 <input type="radio"/>	Spare wipeable inner cover 7150-7227-000 <input type="radio"/>	Spare wipeable inner cover 7150-7228-000 <input type="radio"/>

Side Lying Leg Support

Assists with pelvic and thigh positioning in side-lying

Size 1	Size 2	Size 3
A: 270mm B: 380mm C: 280mm D: 80mm E: 130mm F: 80mm	A: 400mm B: 560mm C: 340mm D: 100mm E: 200mm F: 90mm	A: 400mm B: 610mm C: 400mm D: 110mm E: 150mm F: 70mm
With breathable inner & terry cotton outer 7150-7971-000 <input type="radio"/>	With breathable inner & terry cotton outer 7150-7972-000 <input type="radio"/>	With breathable inner & terry cotton outer 7150-7973-000 <input type="radio"/>
With wipeable incontinence inner & terry cotton outer 7150-7981-000 <input type="radio"/>	With wipeable incontinence inner & terry cotton outer 7150-7982-000 <input type="radio"/>	With wipeable incontinence inner & terry cotton outer 7150-7983-000 <input type="radio"/>
Spare terry cotton cover 7150-7991-000 <input type="radio"/>	Spare terry cotton cover 7150-7992-000 <input type="radio"/>	Spare terry cotton cover 7150-7993-000 <input type="radio"/>

Foot Supports

Provides positioning and skin protection for the feet (pair)

Size 1	Size 2	Size 3	
7150-5351-000 <input type="radio"/>	7150-5352-000 <input type="radio"/>	7150-5353-000 <input type="radio"/>	
A: 21cm (max)	A: 25cm (max)	A: 30cm (max)	

Optional Extras



Stabilising Mesh Set (suitable for foam mattresses only)

Used to stabilise equipment on bed during sleep.
Includes 2x mesh sheets; 1 bottom, 1 top

Cot Bed Size
7150-7011-C00



Single Bed Size
7150-7012-S00





Temperate Topper

Used to allow air to circulate and take moisture away from the body

Cot Bed Size
7150-7051-C00



Single Bed Size
7150-7052-S00




Fitted Sheets

To regulate body temperature and ensure that the client is comfortable

Cot Bed Size
7150-7121-C00



Single Bed Size
7150-7122-S00



SSW Sleep System Carry Bag

7150-CUST-001



Note: for co-sleeping needs, please get in touch with our friendly team.

medifab

Distributed by Medifab Ltd

22 Detroit Drive
Rolleston 7675
New Zealand

0800 543 343

26 Pardoe Street
Devonport TAS 7310
Australia

1300 543 343

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Available in SA from:
Professional Assistance for Living
info@pafl.com.au
08 7092 4555



Script: SimpleStuffWorks_20230725