

Jenx Therapy Bench



THERAPY BENCH

Client Details

| | | |
|-------------------|--|-------|
| Client/User Name: | | Notes |
| Date: | | |
| NDIS/Funding: | | |

Size Options

| Measurements | Size 1 Full Length <input type="radio"/> | Size 1 Half Length <input type="radio"/> | Size 2 Full Length <input type="radio"/> | Size 2 Half Length <input type="radio"/> |
|-------------------------------------------|------------------------------------------|------------------------------------------|------------------------------------------|------------------------------------------|
| Code | 7510-0010-054 | 7510-0015-054 | 7510-0020-054 | 7510-0025-054 |
| Max Load Weight | 85kg | | 85kg | |
| Bench Height from Floor (without cushion) | 165-250mm | | 235-395mm | |
| Bench Top Length | 980mm | 480mm | 980mm | 480mm |
| Bench Top Width | 200mm | | 250mm | |
| Flat Cushion Width | 250mm | | 300mm | |
| Flat Cushion Length | 1000mm | 500mm | 1000mm | 500mm |
| Flat Cushion Depth | 25mm | | 25mm | |


Package includes:

- Height Adjustable Bench
- 1" Comfortable Wipe Down Flat Cushion



Accessories

Half Round Cushion (half length only; order 2 for full length bench)

| Size 1 <input type="radio"/> | Size 2 <input type="radio"/> | |
|-----------------------------------|-----------------------------------|---------------------------------------------------------------------------------------|
| 500mm (L) x 250mm (W) x 125mm (H) | 500mm (L) x 300mm (W) x 150mm (H) |  |
| 7510-0501-054 | 7510-0502-054 | |
| | | |

Pelvic Support

Lateral size 2 support pads with flip away action. Cannot be used with half round cushion.

7510-1900-054



Abduction Block

Available in one size and used to provide simple abduction to ensure the right position for the child.

7510-1750-054



Accessories (cont.)

| | |
|--------------------------------------------|--------------------------------------------|
| <input type="radio"/> Knee Blocks - Size 1 | <input type="radio"/> Knee Blocks - Size 2 |
|--------------------------------------------|--------------------------------------------|

Knee block with fully adjustable knee cups which can be used to control legs and create abduction

| | |
|---------------|---------------|
| 7510-1601-054 | 7510-1602-054 |
|---------------|---------------|



Additional Client Notes

| Notes |
|-------|
| |
| |
| |
| |
| |
| |

medifab®

Distributed by Medifab Ltd

32 Detroit Drive
Rolleston 7675
New Zealand

+64 3 307 9790

solutions@medifab.com

www.medifab.com



DISTRIBUTED BY:

Available in SA from:
Professional Assistance for Living
info@pafl.com.au
08 7092 4555



Script: Jenx_TherapyBench_20181205_Unpriced